VT-1978

Electronic personal scale

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ENGLISH

INSTRUCTION MANUAL ELECTRONIC PERSONAL SCALE

This appliance will help you turn to healthy life-style. Read these instructions carefully and save them for future reference. **The “BIA” technique (bioelectric body impedance analysis) implemented by the appliance allows to calculate fat, water and muscle mass percentage with high accuracy taking into account memorized personal data (height, age, gender).**

WHAT IS EXCESS WIGHT?

It is well known that excess weight may affect the health adversely.

Excess weight may cause fatness, heart problems, progress of atherosclerosis, osteoporosis, although the lack of weight may lead to osteoporosis too (osteoporosis is the systemic skeleton disease. Its features are: bone mass per volume decreases, bone tissue micro architectonics impairment. It issues in bones fragility, the risk of fractures because of minimal trauma or even without one.)

The physical exercises and moderate diet help to get rid of excess weight. Excess weight analyzer is the ideal tool to control excess weight.

It is well known, that in case of small weight muscles and tissues contain more water than fat. It is also common that water conducts electric current better than fat does.

You can measure your body fat percentage in muscles by conducting weak electric current through the body. This helps you control your forms.

Recommendations

- Weigh yourself at the same time every day.
- Percentage of water should be weighed in the evening before supper.
- 2-3 hours should pass after the last eating.
- 10-12 hours should pass after the intensive physical exercises.

Attention!

- **The information about percentage of fat, water and muscles could be obtained only if you stand on the scales barefoot.**
- Percentage of fat depends on different factors and can vary during a day.

Precautions

- The use of excess weight analyzer is

- ABSOLUTELY prohibited for pregnant women.
- The use of excess weight analyzer is ABSOLUTELY prohibited for persons having pacemakers and other medical equipment with accumulators built-in.

FUNCTION BUTTONS

1. ◀ - decrease the number displayed.
2. **ENTER** – enter and confirm the data.
3. ▶ - increase the number displayed.

USAGE

- Open the cover of the batteries compartment and make sure that 2 lithium CR2032 batteries are installed properly.
- Remove isolating gap to connect the batteries.

Selecting units

To select units use the **KG/LB** button situated in the battery compartment.

“**kg**” - weight will be indicated in kilograms and the height in centimeters.;

“**lb**” - weight will be indicated in pounds and the height in feet and inches;

SWITCH ON

- Place the scale on the firm stable surface, if the floor is uneven or covered with carpet, precision will decrease.
- Switch on the unit by slight knock at the case.

Attention!

- If “**oL**” appears on the display, the scale is overloaded (the weight is more than 180 kg.). Immediately get off the scale to prevent damaging of the unit.
- When the “**0000**” symbols are displayed and begin to disappear gradually, the percentage of fat water, muscle mass and bone mass is being measured.
- In case of misuse of the appliance the “**Err2**” (measurements cannot be made, you are not barefoot) will be displayed.
- “**Err %**” – means an error in determination of fat percentage in the muscles.
- “**ErrL**” - percentage of fat in the muscles is so small that measurements cannot be performed.
- “**ErrH**” - percentage of fat in the muscles is so great that measurements cannot be performed.

WEIGHING

- Switch on the unit by slight knock at the case.
- Wait till the **"0.0 kg"** or **"0.0 lb"** symbols, depending on the units selected by the button in the battery compartment, are displayed.
- Stand on the scales, with your feet parallel and evenly distribute the weight.
- Stay motionless while weighing. Flashing symbols are displayed during measuring, after it weight indications will be lighting constantly.
- The **"Err2"** symbols indicates that percentage of water, fat and muscles cannot be determined, (you are not barefoot), in several seconds the scales will switch off automatically.

PERSONAL DATA INPUT AND EXCESS WEIGHT DETERMINATION

1. Switch on the unit by slight knock at the case.
2. Symbols will be displayed and the **"P0"** will start flashing. Enter the number, using the **"◀ / ▶"** buttons. You can input personal data for 10 users P0, P1...P9. Press **"ENTER"** to confirm the user selection.
2. To set the gender press the **"ENTER"** button, the gender indicator will appear:
"♂" – for men, **"♀"** – for women, Press **"◀"** or **"▶"** to change the gender. Confirm your gender selection by pressing the **"ENTER"** button.
3. Press **"ENTER"** once again to set your age data. Input your age by pressing **"◀"** or **"▶"**. The age range is from 10 to 99 years. Press **"ENTER"** to confirm your age input.
4. Press **"ENTER"** once again to set your height data. Input your height by pressing **"◀"** or **"▶"**. The height range is from 100 to 250 cm. Press **"ENTER"** to confirm your height input.
5. If **"0.0 kg"** or **"0.0 lb"** is displayed, stand on the scales barefoot, your feet and heels should be fitted into the circles. Wait until your weight is displayed.
7. Keep staying on the scales. When the **"0000"** symbols are displayed and begin to disappear gradually, the percentage of fat and water, muscles and bone mass in body tissues is being measured.
8. After measurements the following data will be displayed: - percentage of fat **"fat"**, water **"TBW"**, muscles **"MUS"**, bone mass **"BONE"**, as well as daily demand in kilocalo-

rie **"KCAL"**. The results will be repeated three times in this sequence.

9. Estimation of the results is reflected by color indication; – red arrow means excess weight, yellow sector – normal weight, green arrow – underweight.
10. Then weight indications and percentage of fat are recorded into the memory cell of the user, which is displayed by gradually decreasing **"0000"**, **"0.000 - ...0"**, then the scales will switch off.

Note:

- If time range between each personal data input is more than 10 sec, the appliance will switch off automatically.
- Data about the percentage of fat, water, muscles and bone mass can be obtained only if you are standing on the scales barefoot.
- If no any information was provided the appliance would measure fat, water, muscles and bone mass for 25 years old woman 165 cm in height.

Table of percentage of fat

Note: -the data are recommended.

Age	Women Fat %	Men Fat %	Indicators
<=30	4.0-28.0	4.0-22.0	Green
	28.1-33.0	22.1-27.0	Yellow
	33.1-80.0	27.1-80.0	Red
>30	4.0-30.0	4.0-24.0	Green
	30.1-35.0	24.1-29.0	Yellow
	35.1-80.0	29.1-80.0	Red

Changing of fat percentage during a day

Such factors as way of life, physical activity, health status as well as quantity of liquid in the organism influence percentage of fat.

You should take into account the following factors:

- fat percentage can be lower after taking bath or shower.
- fat percentage can be higher after eating or drinking.
- liquid can decrease and thus change fat percentage during or after a disease, and after physical exercises.

ENGLISH

Fat percentage in the organism of a healthy man:

Age	18-39 years		40-59 years		60-79 years	
	Men	Women	Men	Women	Men	Women
Weight deficit	5-7%	14-20%	5-10%	14-22%	5-12%	14-23%
Norm	8-19%	21-32%	11-21%	23-33%	13-24%	24-35%
Excess weight	+20%	+33%	+22%	+34%	+25%	+36%

Source: ISSA (International Sports Science Association)

Water percentage

Water percentage is the amount of water in the body of a man in percentage to total weight.

Water is of primary importance for tissues and organs of a man, about 50-65 % of our organism consist of water. You can figure out percentage of water using a formula, used by dietarians, the result is in percentage from total weight.

Thanks to "BIA" technique used in the scales, you can control percentage of water in your body and adjust this important indicator.

Water norm in the organism for women – 45-60%, for men – 50-65%

What is muscle mass?

Muscle mass includes weight of muscles of supporting-motor system, smooth muscles (heart and digestive system) and water, contained in these muscles.

Muscle mass should be especially controlled by people, suffering from excess weight and obesity. Weight reduction can lead to reduction of fat tissues, but only if muscle mass will be accumulated.

What is bone mass?

Human skeleton is based on bone tissues. They are hard and environmental stable, protecting internal, and together with muscle tissue provide human ability to move.

In the basis of bone mass measurements is weighing of mineral substances in the organism (calcium and others).

According to the investigations, development of muscle tissues strengthens skeleton bones. This fact should be taken into account when composing diet and while physical activity for bones to remain strong.

Standard indications of bone mass for different weight relative age from 20 to 40 years

Women		Men	
Body weight, kgr	Average indicator of bone mass, kg	Body weight, kg	Average indicator of bone mass, kg
Less than 50 kg	1,95 kg	Less than 65 kg	2,66 kg
50 – 75 kg	2,40 kg	65 – 95 kg	3,29 kg
75 kg and more	2,95 kg	More than 95 kg	3,69 kg

Definition of Basal Metabolism Indicator (BMI)

In organism of a human being metabolic processes (breath, blood circulation, digestion, nervous function etc.) always happen, certain energy is consumed, which is expressed in kilocalorie (Kcal). This information allows determining the amount of kilocalorie, which is necessary for metabolic processes. It is minimal energy level, which is necessary for your organism in order to provide normal work of breathing, blood, nervous systems, liver, kidneys and other organs at rest.

Metabolism level increases when you are active. This happens because the supporting-motor system (40% of total organism weight), like a motor, consumes much energy. The greater part of energy is spent in the process of muscles operation; that is why increase of muscle mass allows increasing of Basal Metabolism Indicator.

The higher Basal Metabolism Indicator is, the more calories are spent in the process on muscle mass increasing, and it helps to improve fat level in the organism.

The lower Basal Metabolism Indicator is, the slower fat burning is, and the higher the possibility of obesity with all the consequences that come with it is.

Memory

You can store in memory data for 10 users (P0 ...P9), for each user personal parameters are stored (gender, age, height), as well as two main parameters (weight and fat percentage), which can be viewed for comparison with measuring results.

Access to memory cells

- Switch on the unit by slight knock at the case.
- If **"0.0 kg"** or **"0.0 lb"** are displayed, stand on the scales barefoot, your feet and heels should be fitted into the circles. Wait until your weight is displayed.
- Keep staying on the scales. When the **"0000"** symbols are displayed and begin to disappear gradually, the percentage of fat and water, muscles and bone mass in body tissues is being measured.
- After measurements the following data will be displayed: - percentage of fat **"fat"**, water **"TBW"**, muscles **"MUS"**, bone mass **"BONE"**, as well as daily demand in kilocalorie **"KCAL"**. The results will be repeated three times in this sequence.
- Estimation of the results is reflected by color indication; – red arrow means excess weight, yellow sector – normal weight, green arrow – underweight.
- Then weight indications and percentage of fat are recorded into the memory cell of the user, which is displayed by gradually decreasing **"0000"**, **" . 000 - ...0"**.
- At gradually decreasing of **"0000"**, **" . 000 - ...0"**, you can refer to the necessary memory cell.
- Use **"◀ / ▶"** to select user **"P0, P1...P9"**, and press **"ENTER"** to confirm.
- After it the results of the previous measuring will be displayed (weight and fat percentage), color indicators (red arrow means excess weight, yellow sector – normal weight, green arrow – underweight) will not light up.
- Then gradually decreasing **"0000"**, **" . 000 - ...0"** will appear on the display, and the scales will switch off.

Note:

- Access to memory cells is possible only at the moment of recording of current measuring indications at gradually decreasing **"0000"**, **" . 000 - ...0"** on the display.
- If in one of the memory cells (for instance, **"P5"**) there are no data about weight and fat percentage, it is impossible to recall indications from the cell (for instance **"P6"**).
- If the unit operates unstable, check the batteries and replace them if necessary.

MAINTENANCE

- Take care using the scales. This is the precision tool, avoid temperature differences, high humidity, strikes, dust and direct sun rays.
- Keep away from heaters.
- Do not jump on the scales.
- Wipe the platform with damp cloth.
- Do not use chemicals and abrasive agents to clean the scales.
- Do not let liquids inside the housing of the scales.
- Do not use or keep the scales in humid (more than 80%) cabinets, do not let it to contact with water or other liquids, as such a contacts could affect the measurement accuracy or lead to scale malfunction.
- If you do not use the appliance for a long time, open the battery compartment lid and remove the batteries.
- Keep the appliance in a dry place out of the reach of children.

SPECIFICATIONS

Max.weight: 180 kg

Min.weight: 2 kg

Divisions: 0,1 kg

Battery: 2 batteries 3V CR 2032

Specifications subject to change without prior notice.

Service life of appliance - 5 years

GUARANTEE

Details regarding guarantee conditions can be obtained from the dealer from whom the appliance was purchased. The bill of sale or receipt must be produced when making any claim under the terms of this guarantee.

This product conforms to the EMC-Requirements as laid down by the Council Directive 89/336/EEC and to the Low Voltage Regulation (73/23 EEC)

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