
VT-1983 BK

Personal scale

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PERSONAL SCALE

This model of scale uses bioelectric body impedance analysis (impedance) «BIA», that provides highly accurate calculating of fat and water percentage in human body based on personal data (height, age, sex) saved in the unit memory.

DESCRIPTION

1. Contact pads
2. Data selection sensor buttons «◀ / ▶»
3. Data selection and confirmation sensor button «SET»
4. Battery compartment lid
5. Weight measurement units selection button
6. Display

RECOMMENDATIONS FOR CUSTOMER ATTENTION!

Before using the unit, read this instruction carefully. Keep this instruction for future reference.

Use the unit according to its intended purpose only, as it is stated in this user manual.

Mishandling the unit can lead to its breakage and cause harm to the user or damage to his/her property.

SAFETY MEASURES

- **Using the excess weight analysis function is ABSOLUTELY contraindicated to pregnant women.**
- **Using the excess weight analysis function is ABSOLUTELY contraindicated to people using pacemakers and other medical devices with built-in batteries.**
- Handle your scale with care, as any precision instrument, do not expose it to high or low temperatures, high humidity, never expose it to direct sunlight and do not drop it.
- Use the scale away from heating units.
- Avoid getting of liquids into the scale, this appliance is not waterproof. Do not store or use the scale in a room with high humidity (above 80%), make sure that the scale does not come into contact with water or other liquids as this may cause failure of scale indication or its damage.
- Place the scale on a flat, dry and non-slip surface.
- To increase the indication precision do not use the scale on carpets or rugs.

- Do not step on the scale with wet feet or if the scale surface is damp. You can slip and get injured.
- Be careful while stepping on the scale: step on the scale first with one foot; make sure that the scale is steady and then step on it with both feet.
- Stand on the scale and dispose your weight proportionally. Stand still while weighing. Do not jump on the scale.
- If you are not planning to use the scale for a long time, remove the battery from the battery compartment.
- Strictly follow the polarity when replacing the battery.
- Replace the battery in time.
- Never allow children to use the scale as a toy.
- For children safety reasons do not leave polyethylene bags, used as a packaging, unattended.
- **Attention! Do not allow children to play with polyethylene bags or packaging film. Danger of suffocation!**
- This unit is not intended for usage by children or disabled persons unless they are given all the necessary and understandable instructions by a person who is responsible for their safety on safety measures and information about danger that can be caused by improper usage of the unit.
- Do not repair the scale by yourself.
- Refer to authorized service centers for maintenance.

THE UNIT IS INTENDED FOR HOUSEHOLD USAGE ONLY

IMPORTANT:

Electromagnetic compatibility

Electronic scale can be sensitive to the electromagnetic emission of other units located in close proximity (such as mobile phones, radios, radio model controllers and microwave ovens). In cases of such emission (false or inconsistent data indication on the display), use the unit away from the source of interferences.

BEFORE THE FIRST USE

After unit transportation or storage at low temperature it is necessary to keep it for at least two hours at room temperature before switching on.

ENGLISH

- If there is a protective film on the scale surface, remove it.
- Clean the scale with a soft, slightly damp cloth and then wipe it dry.

BATTERY INSTALLATION AND REPLACEMENT

- Remove the battery compartment lid (4) and insert a "CR2032" battery (supplied with the unit), following the polarity.
- Install battery compartment lid (4) back to its place.

Note: *If there is isolation insert in the battery compartment, remove the battery compartment lid (4), remove the isolation insert and install battery compartment lid (4) back to its place.*

Battery replacement

- When the battery is low, the symbol «Lo» will appear on the display (6).
- Open the battery compartment lid (4), replace the «CR2032» battery with the new one, strictly following the polarity, and close the battery compartment lid (4).
- If you do not use the scale for a long time, remove the battery from the battery compartment.

WHAT IS EXCESS WEIGHT?

It is generally known that the excess weight of the body may adversely affect a person's health. Excess weight can cause obesity, heart problems, atherosclerosis and osteoporosis; short weight may also lead to osteoporosis. (Osteoporosis is a skeleton systemic disease characterized by decrease of a bone mass in relation to its volume as well as malfunction of bone microarchitectonics and leading to increase of brittleness of bones and risk of fracture by minimal trauma or even without it).

Physical exercises and bland diet help to lose excess weight and this scale helps to control your weight.

It is known that at low weight muscles and tissues of the body contain more water than fat, and electric current passes through water better than through fat. Passing weak current through the body the scale quickly and easily calculate the percentage of fat and water in a human body. You can control your weight based on the results of calculation of water and fat percentage.

RECOMMENDATIONS

- Weigh yourself at one and the same time of the day.
- We recommend to measure the percentage of fat and water in the evening before evening meal.
- 2-3 hours should pass after the meal.
- 10-12 hours should pass after intensive physical exercises.
- Fat percentage in the body depends on different factors and can change during the day.
- The results may be inaccurate when calculating fat and water percentage of children under 10 years.
- Professional sportsmen, people suffering from edema and people with nonstandard body proportions can show inaccurate results of fat and water percentage calculation.

Note: *Fat and water percentage data can be calculated only if personal data (sex, age, height) are entered correctly and if you step on the scale with bare feet.*

SWITCHING THE SCALE ON

- Place the scale on a flat solid surface; if the floor is uneven or covered with carpet or rug, the weighing accuracy may decrease.
- Tap on the scale body to switch it on.

MEASUREMENT UNITS SELECTION

- To select measurement units press consequently the button (5) that is located next to the battery compartment lid (4).
- Measurement units are shown on the display (6):
 - «kg» - weight is measured in kilograms, height is measured in centimeters;
 - «lb» - weight in pounds, height in feet and inches;
 - «st:lb» - weight in stones and pounds, height in feet and inches.

WEIGHT MEASUREMENT

- Tap on the scale body to switch it on or touch the «SET» button (3).
- Moving symbols «8...8» will appear on the display (6).
- Wait till the symbols «0.0 kg», «0.0 lb» or «0:0 st:lb» appear depending on the selected measurement units.

- Step on the scale, keep your feet parallel to each other and dispose your weight proportionally.

Stay still while weighing. During weight measurement the weight values will be flashing on the display (6), and after measurement they will be shown constantly.

Attention!

- **If a symbol «oL» appeared on the display (6), it means that the scale is overloaded. Step off the scale immediately to avoid its damage.**
- **The scale will be switched off automatically approximately in 10 seconds after switching on.**

ENTERING OF PERSONAL DATA

- Touch the «SET» button (3), one of the user numbers «P1», «P2»..., «P0» will be flashing on the display (6).
- Select the user number touching the buttons (2) «◀ / ▶». You can enter the data for 10 users marked by symbols «P1», «P2»..., «P0». Confirm your choice touching the «SET» button (3).
- Then the flashing symbol of the sex will be shown on display (6). Select the symbol «♀» for women or «♂» for men, touching the sensor buttons (2) «◀ / ▶». Confirm your choice touching the «SET» button (3).
- Digital values of age will be flashing on the display (6). Set the user's age touching the buttons (2) «◀ / ▶», the age range is from 10 to 80 years with the step 1 year. To save the user's age touch the «SET» button (3).
- Then flashing digital values of height will be shown on the display (6), the digital values of height should be entered in values corresponding to the selected weight measurement units. If you select the weight measurement units «kg», the user's height values are given in centimeters «cm», possible height values are from 100 to 220 cm with the step 1 cm. Set the height values by touching the buttons (2) «◀ / ▶» and then confirm the input data by pressing the «SET» button (3).

CALCULATING WATER AND FAT PERCENTAGE IN A HUMAN BODY

- Touch the «SET» button (3), one of the user numbers «P1», «P2»..., «P0» will be flashing on the display (6).

- Select the necessary user's number «P1», «P2»..., «P0» with the buttons (2) «◀ / ▶».
- The symbols of user's number, sex, age and height will be flashing on the display (6), and then you'll see a moving symbol «8...8».
- Wait till the symbols «0.0 kg», «0.0 lb» or «0:0 st:lb» appear, depending on the selected measurement units.
- Step with your bare feet on the scale so that your feet and heels are on the contact pads (1), stay still during the weight measurement and wait till the weight indications appear.
- Keep standing on the scale. If moving symbols «o...o» appear on the display (6), it means that the scale is calculating the percentage of fat and water in human body.
- Once calculation is finished the percentage of fat/fat and water/TBW and a person's weight will be shown on the display (6). The results of fat and water percentage calculation in human body and the weight of the person will be shown on the display (6) twice.
- Before the scale will be switched off, the percentage of fat and water in human body and the weight of the person are saved into the user's memory cell, this is shown by moving symbols «o...o» on the display (6).

Attention!

- **If a symbol «oL» appeared on the display, it means that the scale is overloaded. Step off the scale immediately to avoid its damage.**
- **If a symbol Err2 appeared on the display, it means that there is an error in fat and water percentage calculation in human body (you may have entered the personal data incorrectly or stepped on the scale with your shoes on).**

MEMORY

This model of scale can keep in its memory 10 users (P1 ... P0), for each user not only his personal data (sex, age, height) are saved but also the results of the last 5 weight measurements (percentage of fat and water in human body and the weight of the person) as well. You can check them to compare with previous results of weight measurement.

Access to memory cells

- During saving the percentage of fat and water indications into the user's memory cell the moving symbols «o...o» are shown

ENGLISH

on the display (6), at this time you can check the results of previous weight measurement.

- During the saving process, when the moving symbols «o...o» are shown on the display, touch the button (2) « \blacktriangleleft / \blacktriangleright » and select one of the previous weight measurements (from 1 to 5) and confirm your selection by touching the «**SET**» button (3).
- The results of the previous weight measurement (weight and percentage of fat and water) will be shown on the display.

Attention:

- **Access to users' memory cells is available only when the current weight measurement results are being saved and the moving symbols «o...o» are shown on the display.**

Note:

- If during the data saving no actions are taken, the scale will be switched off automatically in 10 seconds.
- The data concerning the percentage of fat and water in your body may be obtained only if you stand on the scale with bare feet.
- If no personal data was saved into the user's memory cell, the scale calculates the per-

centage of fat and water in the body of a 25-year old woman with height 165 cm.

Change of fat percentage in human body during a day

The indications of fat percentage in human body are affected by such factors as lifestyle, physical exercises, state of health and water percentage in the body.

You should take into account the following factors:

- the fat percentage indications may be lower if measured after taking a bath or a shower,
- the fat percentage indications may be higher if measured after having meal or drinking,
- percentage of liquid in human body may decrease, affecting fat percentage indications during and after illness and after physical exercises.

Percentage of water in human body

Percentage of water in human body is the ratio of liquid amount in human body given in percents and total weight. Water has primary importance for tissues and organs of the human body, since 50-65% of our body consists of water. With this scale you can control the percentage of water in your body and correct this important characteristic of human body.

Normal percentage of water in the body is 45-60% for women and 50-65% for men

	Age	% of fat			% of water		
		Low	Medium	High	Low	Medium	High
Man	20-39	5.0-7.9	8.0-19.0	19.1-45.0	38.0-54.9	55.0-70.0	70.1-85.0
	40-59	5.0-10.9	11.0-22.0	22.1-45.0	38.0-49.9	50.0-68.0	68.1-85.0
	60-79	5.0-12.9	13.0-25.0	25.1-45.0	38.0-49.9	50.0-65.0	65.1-85.0
Woman	20-39	5.0-20.9	21.0-33.0	33.1-45.0	38.0-47.9	48.0-65.0	65.1-85.0
	40-59	5.0-22.9	23.0-35.0	35.1-45.0	38.0-44.9	45.0-65.0	65.1-85.0
	60-79	5.0-23.9	24.0-36.0	36.1-45.0	38.0-44.9	45.0-65.0	65.1-85.0
Man taking exercises up to 10 hours a week	20-39	5.0-5.9	6.0-12.0	12.1-45.0	38.0-54.9	55.0-71.0	71.1-85.0
	40-59	5.0-5.9	6.0-13.0	13.1-45.0	38.0-49.9	50.0-69.0	69.1-85.0
	60-79	5.0-7.9	8.0-13.0	13.1-45.0	38.0-49.9	50.0-67.0	67.1-85.0
Woman taking exercises up to 10 hours a week	20-39	5.0-13.9	14.0-18.0	18.1-45.0	38.0-47.9	48.0-66.0	66.1-85.0
	40-59	5.0-15.9	16.0-20.0	20.1-45.0	38.0-44.9	45.0-65.0	65.1-85.0
	60-79	5.0-16.9	17.0-20.0	20.1-45.0	38.0-44.9	45.0-65.0	65.1-85.0

Note: *The data given in this table are advisory.*

Temporary connection of the results

Pay attention to the fact that all results of calculations and simple weight measurements will be relevant only at longtime comparison of the results. Short-term changes of the weight within several days are, as a rule, the result of trivial loss of liquid in your body tissues.

Evaluation of the results is given based on changes of total weight, fat and water percentage in body tissues and depends on the duration of the period within which the changes have taken place. Short-term changes of the indications during a day can be significantly different from the changes that took place during a week or a month.

As a rule, short-time changes of indications are connected with percentage of liquid in body tissues, and medium- and long-term changes are caused by changes of fat percentage and muscle mass.

- Thus, short-time weight loss and simultaneous growth or constant fat percentage in body tissues may be caused by decrease of the amount of liquid as the result of physical exercises, visiting sauna or quick weight loss.
- Medium-term increase of weight, decrease of the amount of fat or constant fat percentage in body tissues may be caused by increase of muscle mass.
- If you observe simultaneous decrease of total weight and fat percentage, it means that your diet works, i.e. you lose the fat percentage in your body tissues.
- The best results are achieved if you combine diet with physical exercises. In this case you can observe the medium-term increase of muscle mass in body tissues.

CLEANING AND CARE

- Clean the scale with a soft, slightly damp cloth, and then wipe it dry.
- Never use abrasives or solvents.
- Never immerse the scale into water or other liquids.
- Keep the scale away from children in a dry cool place.

DELIVERY SET

1. Scale – 1 pc.
2. CR2032 battery – 2 pcs. (1pc in the unit)
3. Instruction manual – 1 pc.

TECHNICAL SPECIFICATIONS

Power supply: 1 CR2032 battery, 3 V

Max. capacity: 180 kg.

Min. capacity: 2 kg.

The manufacturer preserves the right to change the design and the specifications of the unit without a preliminary notification.

Unit operating life is 5 years

Guarantee

Details regarding guarantee conditions can be obtained from the dealer from whom the appliance was purchased. The bill of sale or receipt must be produced when making any claim under the terms of this guarantee.

This product conforms to the EMC-Requirements as laid down by the Council Directive 89/336/EEC and to the Low Voltage Regulation (73/23 EEC)

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	Yoshi	Yog' %			Suv %		
		Kam	O'rta	Ko'p	Kam	O'rta	Ko'p
Erkak	20–39	5.0–7.9	8.0–19.0	19.1–45.0	38.0–54.9	55.0–70.0	70.1–85.0
	40–59	5.0–10.9	11.0–22.0	22.1–45.0	38.0–49.9	50.0–68.0	68.1–85.0
	60–79	5.0–12.9	13.0–25.0	25.1–45.0	38.0–49.9	50.0–65.0	65.1–85.0
Ayol							
	20–39	5.0–20.9	21.0–33.0	33.1–45.0	38.0–47.9	48.0–65.0	65.1–85.0
	40–59	5.0–22.9	23.0–35.0	35.1–45.0	38.0–44.9	45.0–65.0	65.1–85.0
60–79	5.0–23.9	24.0–36.0	36.1–45.0	38.0–44.9	45.0–65.0	65.1–85.0	
Bir haftada 10 soatgacha sport bilan shug'ullanadigan erkak							
	20–39	5.0–5.9	6.0–12.0	12.1–45.0	38.0–54.9	55.0–71.0	71.1–85.0
	40–59	5.0–5.9	6.0–13.0	13.1–45.0	38.0–49.9	50.0–69.0	69.1–85.0
60–79	5.0–7.9	8.0–13.0	13.1–45.0	38.0–49.9	50.0–67.0	67.1–85.0	
Bir haftada 10 soatgacha sport bilan shug'ullanadigan ayol							
	20–39	5.0–13.9	14.0–18.0	18.1–45.0	38.0–47.9	48.0–66.0	66.1–85.0
	40–59	5.0–15.9	16.0–20.0	20.1–45.0	38.0–44.9	45.0–65.0	65.1–85.0
60–79	5.0–16.9	17.0–20.0	20.1–45.0	38.0–44.9	45.0–65.0	65.1–85.0	

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GB

A production date of the item is indicated in the serial number on the technical data plate. A serial number is an eleven-unit number, with the first four figures indicating the production date. For example, serial number 0606xxxxxxx means that the item was manufactured in June (the sixth month) 2006.

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