- |  |  |
|--|--|
| <p>Ⓓ <b>Tageslichtlampe</b><br/>Gebrauchsanweisung ..... 2-10</p> <p>ⒼⒷ <b>Daylight lamp</b><br/>Instruction for Use ..... 11-19</p> <p>Ⓕ <b>Lampe à lumière du jour</b><br/>Mode d'emploi ..... 20-28</p> <p>Ⓔ <b>Lâmpada de luz diurna</b><br/>Instrucciones para el uso ..... 29-37</p> | <p>Ⓘ <b>Lampada di luce naturale del giorno</b><br/>Istruzioni per l'uso ..... 38-46</p> <p>Ⓓ <b>Gün ışığı tedavi lambası</b><br/>Kullanma Talimatı ..... 47-55</p> <p>Ⓔ <b>Прибор дневного света</b><br/>Инструкция по применению ..... 56-65</p> <p>Ⓐ <b>Lampa o światle dziennym</b><br/>Instrukcja obsługi ..... 66-74</p> |
|--|--|

## Inhalt

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EG-Richtlinie **2002/96/EC – WEEE** (Waste Electrical and Electronic Equipment).

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## English

### Contents

|    |  |    |
|----|--|----|
| 1  | About the daylight lamp .....              | 12 |
| 2  | Signs and symbols .....                    | 13 |
| 3  | Notes .....                                | 13 |
| 4  | Unit description .....                     | 15 |
| 5  | Initial use.....                           | 15 |
| 6  | Operation .....                            | 16 |
| 7  | Replacing the lamp tubes/<br>starter ..... | 17 |
| 8  | Cleaning and care of the unit ..           | 18 |
| 9  | Disposal .....                             | 19 |
| 10 | Technical specifications.....              | 19 |

### Included in delivery

- Daylight lamp
- These operating instructions

# 1 About the daylight lamp

**Dear valued customer,**

Thank you for choosing one of our products. Our name stands for high-quality, thoroughly tested products for applications in the areas of heat, weight, blood pressure, body temperature, pulse, gentle therapy, massage and air.

With kind regards

Your Beurer team

---

## Why use a daylight lamp?

Haven't we all experienced this?

When the amount of daylight in autumn noticeably decreases or after it has done nothing but rain for days on end, it sometimes doesn't take long before we feel a lack of energy or imbalance creeping up on us, sometimes along with listlessness or even depression. The reason for this can be a light deficiency syndrome associated with the impaired quality of life just described.

Light, and sunlight in particular, is vital to nature and above all also for us humans! Sunlight indirectly controls the production of melatonin, which is released into the blood only in darkness. This hormone tells the body when it is time to sleep. In the months when there is less sunlight, therefore, production of melatonin increases. This makes it harder to get up in the morning, since the body's functions are shut down. If you use the daylight lamp immediately after waking up in the morning—thus as early as possible—you can end the production of melatonin so that your mood takes a positive turn. In addition, light deficiency also hinders production of the happiness hormone serotonin, which decisively affects our well-being. So application of light involves quantitative changes in hormones and neurotransmitters in the brain, which can affect our level of activity,

our feelings and our well-being. Light can help reverse this hormone imbalance.

The simplest and often most effective method of countering light deficiency syndromes is taking long walks in fresh air. Even when the sky is overcast, the light intensity (measured in lux; one lux corresponds to the light of one candle) is still high enough to meet the requirement of at least 2500 lux.

This unit simulates daylight over 10000 lux (15-20 cm) and can be used to treat or prevent light deficiency. A normal electric light, by contrast, is insufficient for building enough serotonin. In a well-lit office, for example, the light intensity is no more than 500 lux. Exercise—outdoors, if possible—also helps counter a bad mood. This is because movement releases endorphins, which cause feelings of happiness. And along with this, we ought above all not to forget: Even autumn and winter have plenty of pleasant sides.

## 2 Signs and symbols

The following symbols appear in these instructions.

**Warning** Warning instruction indicating a risk of injury or damage to health.

**Important** Safety note indicating possible damage to the unit/accessory.

**Note** Note on important information.

## 3 Notes

Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

### Safety notes

#### Warning

- The light of the daylight lamp is only intended for use on the human body.
- Before use, ensure that all packaging materials are removed and that there is no visible damage to the unit or accessories. When in doubt, do not use the unit and contact your dealer or the customer service address provided.
- Also ensure that the daylight lamp has a firm footing.
- The housing of the lamp gets warm during operation. Touching it may result in burns! Always let the lamp cool off before touching it.
- The mandatory minimum distance between the daylight lamp and flammable objects is 1.0 m.
- Connect the unit only to the mains voltage listed on the type plate.
- Do not dip the unit into water and do not use it in wet rooms.
- Do not use it in the vicinity of people who are helpless, children under age 3 or people who are insensitive to heat (people with skin lesions due to illness); do not use after taking painkillers or drugs or drinking alcohol.
- Keep children away from packaging materials (**risk of suffocation**).
- Do not cover up or pack away the unit while it is warm.
- Always remove the mains plug and let the unit cool off before you touch it.
- Do not touch the unit with wet hands while it is plugged in; do not allow any water to be sprayed onto the unit. The unit must be operated only when it is completely dry.
- Protect the unit from strong impacts.
- Do not pull the mains plug from the socket by the cord.
- If the cord or the housing is damaged, contact customer service or the dealer, since special tools are needed to make repairs.
- Disconnection from the mains supply is only guaranteed if the mains plug has been removed from the socket.

## Continued

### General notes

#### Important

- This unit is not intended for commercial or clinical use, but only for individual, private household use!
- If you are taking medication such as pain relievers, antihypertensives or antidepressants, consult your doctor before using the daylight lamp.
- People with retinal diseases as well as diabetics need to be examined by an ophthalmologist before beginning to use the daylight lamp.
- Please do not use if you have an eye disease such as cataracts, glaucoma, diseases of the optic nerve in general and inflammation of the vitreous humour.
- If you have health concerns of any kind, contact your general practitioner!
- Before you use the unit for the first time, remove all packaging materials.
- The lamp tubes and starter are excluded from the warranty.

### Repairs


#### Important

- Only open the unit to replace the tubes and the starter. Please do not attempt to repair the unit yourself. This could result in serious injuries. Failure to observe this regulation shall void the warranty.
- If you need to have the unit repaired, contact customer service or an authorised dealer.

## 4 Unit description

### Overview

- 1 Tubes, 2 x 36 W (inside)
- 2 Transparent cover
- 3 Mains cable
- 4 On/off switch
- 5 Fold-out stand

## 5 Initial use

### Setup

Remove the unit from the packaging. Use the fold-out stand to set it up on a level surface. Select a location that leaves a distance of 20 to 50 cm between the unit and the user. The lamp is most effective within this range.

### Mains connection

Connect the unit only to the mains voltage listed on the type plate.

- Push the plug completely into the mains socket.

#### Note:

- Ensure that there is a socket near the place where the lamp will stand.
- Lay the mains cable in such a way that no one can trip over it.

## 6 Operation

|   |  |
|---|--|
| 1 | <p><b>Switching on the lamp</b></p> <p>► Press the ON/OFF button on the upper edge of the lamp.<br/>The tubes turn on.<br/>ATTENTION!<br/>The brightness of the tubes may be visibly different when switched on. The tubes develop their full luminosity after a couple of minutes.</p>  |
| 2 | <p><b>Enjoying the light</b></p> <p>Sit as close as possible to the lamp, between 20 cm and 50 cm. You can go about your normal activities while using it. You can read, write, make telephone calls, etc.</p> <ul style="list-style-type: none"><li>• Every so often, briefly look directly into the light, since the effect results from the eyes/retina receiving the light.</li><li>• Do not look directly into the light the whole time you are using it, however, since that could lead to overstimulation of the retina.</li><li>• Start with shorter periods of illumination, which you increase over the course of a week.</li></ul> <p><b>Note:</b><br/>Eye aches and headaches may arise after the first times using the lamp. These will go away in later sessions, as the nervous system will become accustomed to the new stimulation.</p> |
| 3 | <p><b>Important instructions</b></p> <p>When using the lamp, maintain the recommended distance of 20 - 50 cm between your face and the lamp.<br/>If the distance is 20 cm, the max. duration of use is 1 hour / day.<br/>If the distance is 50 cm, the max. duration of use is 2 hours / day.<br/>Basically,<br/>the closer you are to the source of light, the less time you should use it.</p>   |
| 4 | <p><b>Enjoying light over longer periods</b></p> <p>Repeat your use of the lamp in the dark seasons for at least 7 successive days, or even longer, depending on your individual needs. If possible, conduct the treatment in the morning hours.</p>   |
| 5 | <p><b>Switching off the lamp</b></p> <p>► Press the ON/OFF button on the upper edge of the lamp.<br/>The tubes turn off. Pull the mains plug from the socket.<br/><b>CAUTION!</b><br/>The lamp is hot after being used. Let the lamp cool off first long enough before you put it away and/or pack it up!</p>  |

## 7 Replacing the lamp tubes/starter

If you determine that your unit has a defective tube, you can replace it yourself. To do so, proceed as follows:


### WARNING:

- First switch off the unit and unplug it.
- Let the unit cool off.
- For safety reasons, do not remove any components or screws other than those specified below.
- Otherwise you are at risk of receiving an electric shock! In addition, failure to observe this regulation shall void the warranty!

1

Lay the unit on a level surface.  
Use a Phillips screwdriver to unscrew the 4 screws on the rear side of the unit.

2


Turn over the unit while holding on to the transparent cover.

3

Carefully remove the cover to replace the tube.


4

Open the tube's two transparent retaining clips by pulling up on them.


5

Carefully pull the tube up over the housing rim (1) and out of the housing, then sideways out of the frame (2).

6

Insert the new tube by following these steps in reverse order.

| Replacing the starter:   |  |  |
|--|--|--|
| 7  | Carry out steps 1 – 3 as described above.  |  |
| 8  | Release the starter by turning it anti-clockwise (1) (< ¼ turn), then pull it up and out (2).  |  |
| 9  | Insert the new starter by following these steps in reverse order.  |  |
| 10   | Before screwing shut the unit, ensure that the retaining clips are again sitting firmly against the tubes and that you have not left any objects (such as packaging material, plastic foil, etc.) in the unit. |  |
|  | Please use only the tubes and starter listed in the technical specifications!  |  |

## 8 Cleaning and care of the unit

The unit should be cleaned occasionally.

### Important

Before you clean the unit, always switch it off, unplug it and let it cool off.

Clean using a damp cloth, to which you can apply a little detergent if necessary.

Do not use abrasive detergents and never immerse the unit in water.

### Important

- Ensure that no water penetrates inside the unit!
- Do not wash the unit in a washing machine!
- Do not touch the unit with wet hands while it is plugged in; do not allow any water to be sprayed onto the unit. The unit must be operated only when it is completely dry.

## 9 Disposal

In the interest of protecting the environment, the unit must not be thrown out with the household waste at the end of its service life.

Dispose of the unit at a suitable local collection or recycling point.

Dispose of the device in accordance with EC Directive

**2002/96/EC – WEEE** (Waste Electrical and Electronic Equipment).

If you have any questions, please contact the local authorities responsible for waste disposal.

## 10 Technical specifications

|  |   |
|--|---|
| Dimensions (WHD)                               | 340 x 480 x 80 mm   |
| Weight   | 2.625 kg  |
| Lighting elements                              | 2 x 36 W tubes (2G 11), Daylight  |
| Nominal data for replacement lighting elements | Philips 2G 11, 36 W,<br>Osram/Dulux 2G 11, 36 W   |
| Nominal data for starter                       | AC 220 – 240 V, 4 to 65 W   |
| Power  | 72 watts  |
| Light intensity                                | 10,000 lux (15-20 cm)   |
| Radiation                                      | Output of radiation beyond the visible spectrum (infrared and UV) is so low that it is harmless to eyes and skin. |
| Mains connection                               | 230 V ~, 50 Hz  |

Subject to technical modifications.

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